



<p>MiniStretcher</p>	<p>Evac+Chair Patient Carrier Seat</p>	<p>Evac+Chair Alarm - Anti-theft Device</p>
<p>PART NO. 1-340</p>  <p>A versatile handling device for lifting patients, bed to bed transfers and for strategic stairway evacuation in conjunction with an Evac+Chair in a life saving emergency situation.</p>	<p>PART NO. 1-325</p>  <p>A portable transfer device for lifting incapacitated persons from a wheelchair into or out of an Evac+Chair, where self transfer is difficult.</p>	<p>PART NO. 1-321 (size - 13cm wide)</p>  <p>A battery operated, key-controlled siren which gives out an audible alarm when the Evac+Chair is detached from its housing.</p>
<p>Evac+Chair Comfy Seat</p>	<p>Evac+Chair Stand</p>	<p>Evac+Chair Wall Mounted Cabinet (300H model only)</p>
<p>PART NO. 1-322C</p>  <p>A cushioned seat which provides additional comfort and support. The Comfy Seat clips to the seat frame of the Evac+Chair and can be purchased as a retrofit.</p>	<p>PART NO. 1-318</p>  <p>A self standing black satin finished support suitable for mounting Evac+Chair, avoiding the need for wall mounting brackets</p>	<p>PART NO. 1-330 Dimensions W 70cm H 115cm D 25cm</p>  <p>A secure epoxy painted metal cupboard for housing Evac+Chair under key-holder management.</p>
<p>Evac+Chair Reception Notice</p>		
<p>PART NO. 1-314REC</p>  <p>A framed notice to visitors to a building reminding them of their responsibilities in an emergency evacuation.</p>		

Evac+Chair International Ltd ParAid House Weston Lane Birmingham B11 3RS UK
0845 230 2253 (local rate - UK only)
+44 (0)121 706 6744 info@evacchair.co.uk Fax +44 (0)121 706 6746
www.evacchair.net

